



# **Thermo Scientific DyNAmo Probe 2-Step qRT-PCR Kit**

## **Technical Manual**

F- 460S	20 cDNA synthesis reactions (20 $\mu$ l each) 100 qPCR reactions (20 $\mu$ l each)
F- 460L	100 cDNA synthesis reactions (20 $\mu$ l each) 500 qPCR reactions (20 $\mu$ l each)



## Table of contents

1. Description.....	4
2. Kit components.....	5
3. Shipping and storage .....	5
4. cDNA synthesis .....	6
4.1 Notes about reaction components .....	6
4.2 Reaction setup .....	9
4.3 cDNA synthesis steps .....	10
5. qPCR.....	11
5.1 Notes about reaction components .....	11
5.2 Reaction setup.....	14
6. Analysis.....	18
6.1 Absolute quantification .....	18
6.2 Relative quantification.....	19
6.3 Melting curve.....	19
7. Troubleshooting.....	20
Appendix I: general molecular biology data.....	22
Product use limitation.....	23
Trademark and patent notices; label licenses.....	23

## 1. Description

Quantitative PCR (qPCR) is a useful technique for the investigation of gene expression, viral load, pathogen detection, and numerous other applications. When analyzing gene expression or viral load, the RNA of interest first needs to be reverse transcribed into cDNA. The following qPCR can either be performed separately (two-step qRT-PCR) or in the same tube with reverse transcription (one-step qRT-PCR).

Thermo Scientific DyNAmo Probe 2-Step qRT-PCR Kit is designed for two-step qRT-PCR of RNA samples from various sources. The kit includes all the necessary reagents for cDNA synthesis and following qPCR. Either total RNA, messenger RNA, viral RNA or in vitro transcribed RNA can be used as a template for reverse transcription. The kit includes both random primers and oligo(dT) primers. The user can choose either of these or alternatively use gene specific primers.

The reverse transcriptase in the DyNAmo™ Probe 2-Step qRT-PCR Kit is M-MuLV RNase H<sup>+</sup>, which provides higher sensitivity to qPCR than RNase H<sup>-</sup> reverse transcriptases. The RNase H activity in the RT enzyme facilitates annealing of PCR primers to the cDNA by degrading the RNA template before the PCR step.

The performance of the qPCR step is based on a hot start *Thermus brockianus* (*Tbr*) DNA polymerase. *Tbr* DNA polymerase is chemically engineered to be inactive at room temperature. The inactivation prevents the extension of nonspecifically bound primers during reaction setup and therefore increases PCR specificity. The reaction setup can be performed at room temperature. The initial denaturation step in the PCR protocol reactivates the polymerase (hot start).

The reaction chemistry of DyNAmo Probe 2-Step qRT-PCR Kit is applicable to most block-based and capillary-based real-time PCR instruments, including those from Applied Biosystems, Roche, Bio-Rad Laboratories, Corbett Research and Stratagene. For SYBR® Green based detection, we recommend Thermo Scientific DyNAmo SYBR Green qPCR Kits.

## 2. Kit components

cDNA synthesis	F-460S	F-460L
M-MuLV RNase H <sup>+</sup> reverse transcriptase (includes RNase inhibitor)	1 x 40 µl (sufficient for 20 RT reactions of 20 µl)	1 x 200 µl (sufficient for 100 RT reactions of 20 µl)
2x RT buffer (includes dNTP mix and 10 mM MgCl <sub>2</sub> *)	1 x 200 µl	1 x 1 ml
Random hexamers (300 ng/µl)	1 x 20 µl	1 x 100 µl
Oligo(dT) <sub>15</sub> primer (100 ng/µl)	1 x 20 µl	1 x 100 µl
qPCR	F-460S	F-460L
2x master mix (hot start <i>Tbr</i> DNA polymerase, optimized PCR buffer, MgCl <sub>2</sub> , dNTP mix including dUTP)	1 x 1 ml (sufficient for 100 reactions of 20 µl or 40 reactions of 50 µl)	5 x 1 ml (sufficient for 500 reactions of 20 µl or 200 reactions of 50 µl)
50x ROX passive reference dye	1 x 50 µl	1 x 250 µl

\* Provides 5 mM MgCl<sub>2</sub> in 1x reaction concentration

Material safety data sheet (MSDS) is available at [www.thermoscientific.com/fzmsds](http://www.thermoscientific.com/fzmsds).

## 3. Shipping and storage

The DyNAmo Probe 2-Step qRT-PCR Kit is shipped on gel ice. Upon arrival, store all kit components at -20°C. When using the 2x master mix, the leftover thawed mix can be refrozen and stored at -20°C without affecting the performance of the kit.

## 4. cDNA synthesis

### 4.1 Notes about reaction components

Table 1. General recommendations.

Categories	Comments
Reaction volume	20 $\mu$ l
Template amount	Max 1 $\mu$ g of RNA.
Priming options	Random hexamers, oligo(dT) or gene-specific primers.

#### RNA template

Total RNA, mRNA, viral RNA or *in vitro* transcribed RNA can be used as a template. Successful cDNA synthesis is dependent on the integrity and purity of the template RNA. RNA preparation should be free of any DNA or RNase contamination. The purity of RNA can be determined by measuring the ratio of  $A_{260}/A_{280}$ . The optimal ratio is 1.8–2.0.

RNA isolation should be performed under RNase-free conditions. Furthermore, any contamination with RNases from other potential sources like glassware, plasticware and reagent solutions has to be avoided. This can be done by wearing gloves and using sterile tubes and pipet tips. Water used for the reactions should also be RNase free, but not DEPC treated as traces of DEPC can inhibit PCR.

DNA contamination can be removed from the RNA sample by treating the sample with RNase free DNase I. This should be done especially if primers or probe(s) for the qPCR step cannot be designed in exon-exon boundaries or in separate exons. See information about qPCR primer design in Section 5.1.

The maximum amount of template RNA for DyNAmo Probe 2-Step qRT-PCR Kit is 1  $\mu$ g. This amount includes all RNA present in the sample, for example carrier RNA and other possible RNAs in addition to the target RNA.

#### M-MuLV RNase H<sup>+</sup> reverse transcriptase

M-MuLV RNase H<sup>+</sup> RT provides good sensitivity in qRT-PCR applications, where amplicons are usually around 100 bp in length. Also, with M-MuLV RNase H<sup>+</sup> there is no need to perform separate RNase H treatment after cDNA synthesis, as the RNase H activity in the enzyme degrades RNA in the RNA-cDNA hybrid.

#### RNase inhibitor

The RNase inhibitor included in the mix with the reverse transcriptase inhibits contaminating RNases present in the RNA sample. It does not affect the RNase H activity in the M-MuLV reverse transcriptase.

## RT primers

Specific primers, random hexamers or oligo(dT) primers can be used for the RT step. Using specific primers can help to decrease background, whereas random hexamers and oligo(dT) primers are useful if several different amplicons need to be analyzed from a small amount of starting material.

A good starting point is to use random hexamers for cDNA synthesis. Random hexamers transcribe all RNA (mRNA, rRNA, tRNA and *in vitro* transcribed RNA) producing cDNA that covers the whole transcript. The recommended amount of random hexamers per 20 µl RT reaction is 300 ng (can be optimized between 200–400 ng if necessary).

Oligo(dT) primers can be used to transcribe poly(A)+ RNAs. These include eukaryotic mRNAs and retroviruses with poly(A)+ tails. Several different mRNAs are transcribed allowing subsequent qPCR detection of different targets from the same cDNA synthesis reaction. The recommended amount of oligo(dT) primers per 20 µl RT reaction is 100 ng (can be increased up to 1 µg if necessary). If the amplicon is located at the 5' end of the transcript, using random hexamers is recommended.

Gene-specific primers are used to transcribe only the particular RNA of interest, in contrast to oligo(dT)/random primers that transcribe all mRNAs/RNAs in the sample. The recommended amount of specific primer per 20 µl RT reaction is 10 pmol (can be optimized between 5–20 pmol if necessary).

## Minus RT control

A minus RT control should be included in all qRT-PCR experiments to test for DNA contamination (such as genomic DNA or PCR product from a previous run). Such a control reaction contains all the reaction components except for the reverse transcriptase. RT reaction should not occur in this control, so if PCR amplification is seen, it is most likely derived from contaminating DNA.

## RT efficiency

The cDNA synthesis step is very critical in qRT-PCR. The efficiency of reverse transcription varies and can be low in some cases. The expression level of the target RNA molecule and the efficiency of the RT reaction must therefore be considered when determining the appropriate amount of starting template for subsequent PCR steps. The volume of cDNA template should not exceed 10 % of the qPCR reaction volume, as elevated volumes of template may reduce the efficiency of the PCR amplification. A dilution series of the template can be done to optimize the volume of the starting material used.

## Standards

The absolute amount of the target nucleic acid (expressed as a copy number or concentration) is determined by comparison of Cq values to external standards containing a known amount of nucleic acid. (Cq = quantification cycle, the fractional PCR cycle at which the target is quantified in a given sample. The level of Cq is set manually or calculated automatically.) The external standards should contain the same or nearly the same sequence as the template of interest. It is especially important that the primer binding sites are identical to ensure equivalent amplification efficiencies of both standard and target molecules. In capillary instruments, Cp (crossing point, the cycle number at detection threshold) replaces the term Cq, which is used with block-based instruments.

Using RNA molecules as standards for RNA quantification is recommended. The use of RNA standards takes the variable efficiency of the reverse transcription into account. RNA standards can be generated for example by cloning the cDNA of interest to a vector containing RNA polymerase promoter, e.g. T7 or Sp6. From the vector the insert can be *in vitro* transcribed to obtain the final RNA standard with identical sequence to the target amplicon. The vector must then be degraded with RNase-free DNase, and the concentration of the RNA standard determined spectrophotometrically. Alternatively a defined RNA preparation, e.g. from a cell line or a virus, with known concentration can be used as an RNA standard.

## Reference genes

When studying gene expression, the quantity of the target gene transcript needs to be normalized against variation in the sample quality and quantity between samples. To ensure identical starting conditions, the relative expression data have to be normalized with respect to at least one variable, such as sample size, total amount of RNA, or reference gene(s), for example. A gene used as a reference should have a constant expression level that is independent of the variation in the state of the sample tissue. Examples of commonly used reference genes are beta actin, GAPDH and 18S rRNA. A problem is that, even with housekeeping genes, the expression usually varies to some extent. That is why several reference genes are usually required, and their expression needs to be checked for each experiment.

The amplification efficiency of a reference gene should be the same as the amplification efficiency of the target gene. If this is not the case, the results have to be corrected for the efficiency.

Since RNA quantification involves a number of variables, and each experiment is inherently different, careful experimental design is very important. Useful information and guidelines for experimental design, normalization, RNA standards, etc. can be found in the following review articles:

Bustin S.A. (2000) *Journal of Molecular Endocrinology* 25, 169–193

Bustin S.A. (2002) *Journal of Molecular Endocrinology* 29, 23–39.

## 4.2 Reaction setup

- Perform the reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- All plasticware should be RNase-free.
- Use gloves to prevent RNase contamination.
- Pipette with sterile filter tips.
- Minimize pipetting errors by using calibrated pipettes and by preparing premixes to avoid pipetting very small volumes.
- Pipet all components on ice.
- Reaction tubes should be centrifuged before starting the incubations to force the solution to the bottom of the tubes and to remove any bubbles.

### Protocol

1. Program the cycler as outlined in Table 3.
2. Thaw template RNA, 2x RT buffer and primers. Mix the individual solutions to ensure homogeneity.
3. Prepare a cDNA synthesis premix by mixing 2x RT buffer, primers, RNase free H<sub>2</sub>O and reverse transcriptase (see Table 2). Mix thoroughly to ensure homogeneity. Dispense appropriate volumes into reaction tubes.
4. Add template RNA to the reaction tubes.
5. Place the tubes in the thermal cycler and start the program.

Table 2. Reaction setup.

Components	Stock	20 µl reaction	Comments
RT buffer	2x	10 µl	RT buffer includes dNTPs and MgCl <sub>2</sub>
Random hexamer primer set	300 ng/µl	1 µl	Alternatively oligo(dT) primer or a specific primer can be used. See Section 4.1 ('RT primers').
Template RNA		X µl	Max 1 µg.
M-MuLV RNase H <sup>+</sup> reverse transcriptase		2 µl	Includes RNase inhibitor.
RNase free H <sub>2</sub> O		add to 20 µl	

Table 3. Cycling protocol for reverse transcription.

Step	Purpose	Temperature	Time	Comments
1	Primer extension	25°C	10 min	This step is not necessary if gene specific primers are used.
2	cDNA synthesis	37°C	30 min	Most targets can be synthesized at 37°C. The temperature can be varied between 37–48°C if necessary. Incubation time can be extended up to 60 min if needed for long or rare transcripts.
3	Reaction termination	85°C	5 min	Inactivation of M-MuLV prevents it from inhibiting qPCR reaction.
4	Cooling the sample	4°C	Hold	Optional

### 4.3 cDNA synthesis steps

#### Pre-denaturation (optional)

A separate RNA denaturation step is generally not required, but it can be performed before cDNA synthesis if the template RNA has a high degree of secondary structure. The denaturation step, at 65°C for 5 min, should be performed before adding 2x RT buffer and reverse transcriptase to the reaction mix.

#### Primer extension

The incubation of at 25°C for 10 min extends random primers or oligo(dT) primers before the actual cDNA synthesis. Without the incubation at 25°C the primers may dissociate from the template when the temperature is raised to the cDNA synthesis temperature. This preliminary extension step is not necessary for gene-specific primers.

#### cDNA synthesis

Incubation at 37°C will work for most templates, but it can be optimized between 37°C and 48°C if necessary. Raising the temperature can be helpful if the template has strong secondary structures. Higher temperature can also improve specificity if gene-specific primers are used. Incubation time of 30 min is sufficient in most cases. If the target is located near the 5' end of a long transcript and oligo(dT) priming is used, or the target is rare, cDNA synthesis time can be extended up to 60 min.

#### Reaction termination

The termination step at 85°C inactivates the M-MuLV reverse transcriptase, thus preventing it from inhibiting the qPCR reaction.

## 5. qPCR

### 5.1 Notes about reaction components

Table 4. General recommendations for qPCR.

Categories	Comments
Consumables	Follow the recommendations of the PCR instrument manufacturer.
Reaction volume	20–50 $\mu$ l
Amplicon size	50–250 bp
Template amount	The volume of cDNA template should not exceed 10 % of the qPCR reaction volume.
Primer design	Use primers with matched $T_m$ . Avoid inter-primer and intra-primer complementary sequences. We recommend calculating $T_m$ by the nearest-neighbor method as described by Breslauer <i>et al.</i> (1986) <i>Proc. Nat. Acad. Sci.</i> 83, 3746-50. Instructions for $T_m$ calculation and a link to a calculator using a modified nearest-neighbor method can be found on the Thermo Scientific website ( <a href="http://www.thermoscientific.com/pcrwebtools">www.thermoscientific.com/pcrwebtools</a> ).

#### Probe-based detection chemistries

Many qPCR chemistries based on the use of labeled probes have been developed. Usually the probe is labeled with a fluorophore, and fluorescence of the probe is changed as a consequence of its annealing to the target DNA.

Hydrolysis probe chemistry (TaqMan<sup>®</sup>, Double Dye, etc.) is the most widely used probe-based chemistry in real-time PCR. A hydrolysis probe consists of a target specific sequence, which is usually around 20 bp long. The probe has a fluorescent reporter molecule (fluorophore) in one end and a quencher in the other end of the probe. The quencher receives the energy from the fluorophore and quenches the fluorescence. During the PCR protocol the probe hybridizes to its complementary sequence in the target and one of the PCR primers anneals to the same strand close upstream from the probe. When the polymerase extends the primer it encounters the probe, hydrolyses it from the 5' end, and thus cleaves the reporter from the probe. When the reporter is cleaved it is no more quenched and the increase in the fluorescence can be measured with the real-time PCR instrument.

#### cDNA template

If the cDNA synthesis reaction will not be used for qPCR immediately, it can be stored at -20°C. Also, if only part of the reaction volume is needed for qPCR, store the remainder at -20°C.

The volume of cDNA template should not exceed 10 % of the qPCR reaction volume, as elevated volumes of template may reduce the efficiency of the PCR amplification. Excess salt and random primers in the cDNA synthesis reaction can inhibit the DNA polymerase. A dilution series of the template can be made to optimize the amount of the starting material used.

## Hot start *Tbr* DNA polymerase

The hot start *Tbr* DNA polymerase is a chemically reversibly inactivated enzyme. The inactivation prevents the extension of nonspecifically bound primers during reaction setup and the first heating cycle, and therefore increases PCR specificity. The initial denaturation step in the PCR protocol reactivates the polymerase (hot start). Due to the hot start polymerase, the reaction setup can be performed at room temperature. The hot start *Tbr* DNA polymerase has 5'→3' exonuclease activity, which is required for hydrolysis probe chemistries, e.g. for TaqMan® chemistry.

## Primers and probe(s) for qPCR step

Careful primer and probe design is particularly important to minimize nonspecific primer annealing and primer-dimer formation. Standard precautions must be taken during primer design to avoid primer-dimer or hairpin loop formation. Many software tools for designing PCR primers and probes simultaneously are available.

The optimal concentration for primers is usually between 0.05 and 1 µM and for probe between 0.05 and 1 µM, but the optimum depends on the chemistry and other assay variables. Requirements for probe design and probe concentration depend on the chemistry used. E.g. for TaqMan chemistry the recommended starting concentration for primers is 0.5 µM and for probe 0.25 µM.

PCR primers in qRT-PCR experiments should be designed to anneal to sequences in two exons on opposite sides of an intron. A long intron inhibits the amplification of the genomic target. Alternatively, primers or probe(s) can be designed to anneal to the exon-exon boundary of the mRNA. With such an assay design, the priming of genomic target is highly inefficient.

## ROX™ passive reference dye

For most real-time instruments ROX™ passive reference dye is not required, but on some instruments it is used to normalize for non-PCR related fluorescence signal variation. Passive reference dye does not take part in the PCR reaction and its fluorescence remains constant during the PCR reaction. The amount of the ROX passive reference dye needed can vary depending on the type of the excitation. The amount of ROX dye needed with real-time cyclers which use argon laser as the excitation light source or which have excitation filters that are not optimal for ROX dye may be greater than with instruments that excite efficiently near 585 nm.

The ROX dye is provided as a 50x solution dissolved in a buffer that is compatible with the qPCR reaction buffer. The optimal ROX dye concentration is usually 0.3–1x (see Table 5 for instrument-specific recommendations). Note that the use of ROX passive reference dye may not be possible with some fluorescent dyes.

Table 5. ROX concentration.

Real-time PCR instrument	Recommended ROX concentration
Applied Biosystems StepOne™ Real-Time PCR System	1x
Applied Biosystems 7000, 7300, 7700 Real-Time PCR Systems	1x
Applied Biosystems 7900HT Real-Time PCR System	1x
Applied Biosystems 7500 Real-Time PCR System	0.3x
Agilent Mx3000P® QPCR System	0.3x (optional)
Agilent Mx3005P® QPCR System	0.3x (optional)
Agilent Mx4000® QPCR System	0.3x (optional)

### UNG (UDG) treatment

Due to the high sensitivity of qPCR, even minute amounts of contaminating DNA can lead to false positive results. If dUTP is used in all qPCR reactions, the carry-over contamination from previous PCR runs can be prevented by treating the reaction samples with UNG before PCR. UNG (uracil-N-glycosylase) digests dU-containing DNA, and the digested DNA cannot act as a template in qPCR (Longo M.C. *et al.* (1990) *Gene* 93: 125–28). UNG is inactivated during the first denaturation step in PCR. The UNG treatment step (50°C for 2 min) has no negative effect on qPCR performance because the hot-start *Tbr* DNA polymerase is not reactivated at 50°C. All Thermo Scientific DyNAmo qPCR Kits contain dUTP and therefore UNG treatment can be used.

To minimize contamination risk in general, tubes containing reaction products should not be opened or analyzed by gel electrophoresis in the same laboratory area that is used to set up reactions.

## 5.2 Reaction setup

- Perform the reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- As the hot-start *Tbr* DNA polymerase is inactive during PCR setup, it is not necessary to do the setup on ice.
- Pipette with sterile filter tips.
- Minimize the exposure to light after adding ROX passive reference dye and/or probe to the 2x master mix.
- Minimize pipetting errors by using calibrated pipettes and by preparing premixes to avoid pipetting very small volumes.
- Use optically clear caps or sealers to achieve maximum signal.
- Use a cap sealing tool or firm finger pressure to close caps properly, or use a film sealer.
- Avoid touching the optical surface of the cap or sealing film without gloves, as fingerprints may interfere with fluorescence measurements.
- Plates or strips should be centrifuged before starting the cycling program to force the solution to the bottom of the tubes and to remove any bubbles.
- Use molecular biology grade H<sub>2</sub>O.

### 5.2.1 General protocol for all instruments

If you are using an Applied Biosystems real-time PCR instrument, see Section 5.2.2.

#### Reaction setup

1. Program the cycler as outlined in Table 7.
2. Thaw the template cDNA, primers, probe and master mix (and the ROX passive reference dye, if necessary). Mix the individual solutions to ensure homogeneity. This is especially important for the master mix.
3. Prepare a PCR premix by mixing the master mix, primers, probe (ROX if used,) and H<sub>2</sub>O. Mix the PCR premix thoroughly to ensure homogeneity. Dispense appropriate volumes into strip tubes or plate wells.
4. Add template cDNA to the strip tubes or plate wells containing the PCR premix. The volume of the cDNA added (from the RT reaction) as the template should not exceed 10 % of the final PCR volume.
5. Seal the strips or plate with appropriate sealer, place them in the thermal cycler and start the cycling program.

Table 6. Reaction setup for Hydrolysis probes (TaqMan, Double Dye, etc.).

Components (In order of addition)	50 $\mu$ l reaction	20 $\mu$ l reaction	Final concentration	Comments
2x Master mix	25 $\mu$ l	10 $\mu$ l	1x	Mix thoroughly. Avoid air bubble formation.
Primer mix (in H <sub>2</sub> O)	X $\mu$ l	X $\mu$ l	0.5 $\mu$ M fwd 0.5 $\mu$ M rev	Titrate from 0.05 to 1 $\mu$ M if necessary.
Probe	X $\mu$ l	X $\mu$ l	0.25 $\mu$ M (TaqMan probe)	Titrate from 0.05 to 0.5 $\mu$ M if necessary.
50x ROX reference dye	(0.03–1 $\mu$ l)	(0.012–0.4 $\mu$ l)	0.03–1x	Optional. See Section 5.1 and 5.2.2.
Template cDNA	X $\mu$ l	X $\mu$ l		Do not exceed 10 % of the final reaction volume. A dilution series of the cDNA synthesis reaction can be made to optimize the amount.
H <sub>2</sub> O	add to 50 $\mu$ l	add to 20 $\mu$ l		Add water to fill up to the final reaction volume.

For different volumes, adjust all components proportionally.

## Cycling protocol

Table 7. Cycling protocol for Hydrolysis probes (TaqMan, Double Dye, etc.).

Step	Purpose	Temp.	Time	Comments
	UNG incubation			Optional, see below.
1	Initial denaturation	95°C	15 min	This step is needed to activate the hot start <i>Tbr</i> DNA polymerase and to denature the template cDNA.
2	Denaturation	95°C	15 s	
3	Annealing + extension	60°C	60 s	Temperature and time can be adjusted, but it is important that the polymerase and exonuclease activities are both functional at this step and the probe is hybridized to the target sequence.
4	Data acquisition			Fluorescence data collection
5	Number of cycles	35–45 cycles, steps 2–4		

### **UNG incubation (optional)**

If UNG enzyme is used, incubate 2 min at 50°C. This step does not negatively affect qPCR performance, because the hot-start DNA polymerase is not active at 50°C. If heat-labile UNG is used, decrease the incubation temperature and increase time in accordance with the manufacturers' instructions.

### **Initial denaturation / reactivation**

Initial denaturation at 95°C for 15 minutes is needed to ensure complete reactivation of the hot-start DNA polymerase and denaturation of the template.

### **Denaturation**

Denaturation at 95°C for 10 s is sufficient in most cases.

### **Annealing/extension**

For most amplicons, a combined annealing and extension for 60 seconds at 60°C works well if the primers are designed to anneal efficiently at 60°C ( $T_m$  about 65°C). An annealing temperature of 60°C has proven successful with a wide range of primer pairs.

These guidelines are based on  $T_m$  values (50 mM salt and 0.5  $\mu$ M primer) calculated by the nearest-neighbor method as described by Breslauer *et al.* (1986) *Proc. Nat. Acad. Sci.* 83: 3746–50. Instructions for  $T_m$  calculation and a link to a calculator using a modified nearest-neighbor method can be found on the Thermo Scientific website ([www.thermoscientific.com/pcrwebtools](http://www.thermoscientific.com/pcrwebtools)). Different software may give different  $T_m$  values.

If needed, the annealing temperature can be optimized by performing additional runs, varying the annealing temperature in each by 2°C. A temperature gradient feature on the thermocycler can also be used, if available.

### **Number of cycles**

For most applications, 40 cycles of amplification should be sufficient even when the template is present at a very low copy number.

## 5.2.2 Protocol for Applied Biosystems real-time PCR instruments requiring ROX

### Addition of ROX™ passive reference dye


#### ABI 7000, 7300, 7700, 7900 and StepOne™: 1x ROX final concentration

1. Thaw and carefully mix 50x ROX and 2x Master mix tubes.
2. Add 40 µl of 50x ROX to each 1 ml 2x Master mix tube.
3. Mix again carefully.
4. Store at -20°C.



#### ABI 7500: 0.3x ROX final concentration

1. Thaw and carefully mix 50x ROX and 2x Master mix tubes.
2. Add 12 µl of 50x ROX to each 1 ml 2x Master mix tube.
3. Mix again carefully.
4. Store at -20°C.



### Reaction setup for all Applied Biosystems models:

1. Program the cyclers as outlined in Table 9.
2. Thaw template DNA, primers, probe(s) and 2x Master mix (where ROX passive reference dye has been added). Mix the individual solutions to assure homogeneity. This is especially important for the Master mix.
3. Prepare a PCR premix by mixing 2x Master mix, primers, probe(s), and H<sub>2</sub>O. Mix the PCR premix thoroughly to ensure homogeneity. Dispense appropriate volumes into strip tubes or plate wells. Use reverse pipeting technique to avoid bubbles.
4. Add template DNA (<200 ng/20 µl reaction) to the strip tubes or plate wells containing the PCR premix. For two-step qRT-PCR, the volume of the cDNA added (from the RT reaction) should not exceed 10 % of the final PCR volume.
5. Seal the strips or plate with appropriate sealer, place them in the thermal cycler and start the cycling program.

Table 8. Reaction setup for Applied Biosystems real-time PCR instruments.

Components (In order of addition)	50 $\mu$ l reaction	20 $\mu$ l reaction	Final concentration	Comments
2x Master mix with ROX added (see instructions above)	25 $\mu$ l	10 $\mu$ l	1x	Mix thoroughly. Avoid air bubble formation.
Primer mix (in H <sub>2</sub> O)	X $\mu$ l	X $\mu$ l	0.3 $\mu$ M fwd 0.3 $\mu$ M rev	Titrate from 0.05 to 1 $\mu$ M, if necessary
Probe	X $\mu$ l	X $\mu$ l	0.25 $\mu$ M (TaqMan <sup>®</sup> probe)	Titrate from 0.05 to 0.5 $\mu$ M, if necessary.
Template DNA (in H <sub>2</sub> O)	X $\mu$ l	X $\mu$ l		Do not exceed 10 ng/ $\mu$ l in the final reaction.
H <sub>2</sub> O	add to 50 $\mu$ l	add to 20 $\mu$ l		

Table 9. Cycling protocol for Applied Biosystems real-time PCR instruments.

Step	Temp.	Time	Cycles
Initial denaturation	95°C	15 min	1
Denaturation	95°C	15 s	40 cycles
Annealing/extension	60°C	60 s	

## 6. Analysis

### 6.1 Absolute quantification

For RNA quantification, the use of standards is recommended (see 'Standards' in Section 4.1). Absolute quantification is performed by plotting samples of unknown concentration on a standard curve generated from a dilution series of template RNA of known concentration. Typically, the standard curve is a plot of the quantification cycle (C<sub>q</sub>) against the logarithm of the amount of RNA. A linear regression analysis of the standard plot is used to calculate the amount of RNA in unknown samples. The slope of the equation is related to the efficiency of the PCR reaction. The PCR efficiency should be the same for standards and samples for quantification to be accurate. The PCR efficiency of the samples can be determined by doing a dilution series of these samples.

For a graph where Cq is on the y axis and log(DNA copy #) on the x axis:

$$\text{PCR efficiency} = \left(10^{\frac{-1}{\text{slope}}}\right) - 1 \times 100 \%$$

A slope of -3.322 corresponds to 100 % efficiency.

For a graph where log(DNA copy#) is on the y axis and Cq on the x axis:

$$\text{PCR efficiency} = \left(10^{-1 \times \text{slope}}\right) - 1 \times 100 \%$$

A slope of -0.301 corresponds to 100 % efficiency.

## 6.2 Relative quantification

Relative quantification is used to determine the ratio between the quantity of a target molecule in a sample and in the calibrator (healthy tissue or untreated cells, for example). The most common application of this method is the analysis of gene expression, such as comparisons of gene expression levels in different samples, for example. The target molecule quantity is usually normalized with a reference gene (see 'Reference genes' in Section 4.1).

If the amplification efficiency of a reference gene is the same as that of the target gene, the comparative  $\Delta\Delta\text{Cq}$  method can be used for relative quantification. Both the sample and the calibrator data are first normalized against variation in sample quality and quantity. Normalized ( $\Delta\text{Cq}$ ) values are calculated by the following equations:

$$\Delta\text{Cq}(\text{sample}) = \text{Cq}(\text{target}) - \text{Cq}(\text{reference})$$

$$\Delta\text{Cq}(\text{calibrator}) = \text{Cq}(\text{target}) - \text{Cq}(\text{reference})$$

The  $\Delta\Delta\text{Cq}$  value is then determined using the following formula:

$$\Delta\Delta\text{Cq} = \Delta\text{Cq}(\text{sample}) - \Delta\text{Cq}(\text{calibrator})$$

The expression of the target gene normalized to the reference gene and relative to the calibrator =  $2^{-\Delta\Delta\text{Cq}}$

If the amplification efficiency of a reference gene is not the same as that of the target gene, a method should be used that takes this into account (Pfaffl MW., (2001) *Nucleic Acids Res.* 29: e45).

## 6.3 Melting curve

Melting curve analysis is possible only with some probe chemistries. It cannot be performed e.g. with TaqMan chemistry. Melting curve analysis is typically included in the analysis software of real-time fluorescence detection instruments. Data interpretation depends on the chemistry and the application.

## 7. Troubleshooting

Possible causes	Comments and suggestions
<b>No increase in fluorescence signal</b>	
Error in cycler setup	<ul style="list-style-type: none"> <li>Make sure that the instrument settings are correct for the experiment.</li> </ul>
Missing components (e.g. primers, probe or template) or pipetting error	<ul style="list-style-type: none"> <li>Check the assembly of the reactions.</li> <li>Check the concentrations and storage conditions of the reagents.</li> </ul>
RNA degraded or poor quality	<ul style="list-style-type: none"> <li>Check the concentration, integrity, purity and storage conditions of the RNA template. Make new RNA dilutions from the stock if necessary.</li> </ul>
Incorrect temperature in cDNA synthesis reactions	<ul style="list-style-type: none"> <li>The recommended temperature in cDNA synthesis step is 37°C. It can be optimized between 37–48°C if necessary.</li> </ul>
Missing essential step in the cycler protocols	<ul style="list-style-type: none"> <li>Check the cycler protocols for cDNA synthesis and qPCR steps.</li> </ul>
RT-PCR product too long	<ul style="list-style-type: none"> <li>The length of the amplicon should be between 50 and 250 bp. The optimal length is 100–150 bp.</li> </ul>
qPCR primer/probe design or concentration not optimal	<ul style="list-style-type: none"> <li>Check primer/probe design and concentration. See Section 5.1.</li> </ul>
Sample not configured properly in the cycler software	<ul style="list-style-type: none"> <li>Check the plate configuration fed into the cycler software.</li> </ul>
<b>Late increase in fluorescence signal</b>	
Error in cycler setup	<ul style="list-style-type: none"> <li>Make sure that the instrument settings are correct for the experiment.</li> </ul>
Missing components or pipetting error	<ul style="list-style-type: none"> <li>Check the assembly of the reactions.</li> <li>Check the concentrations and storage conditions of the reagents.</li> </ul>
RNA template amount too low	<ul style="list-style-type: none"> <li>Check the calculation of the template stock concentration; increase the amount of RNA template (max 1 µg).</li> </ul>
RNA degraded or poor quality	<ul style="list-style-type: none"> <li>Check the concentration, integrity, purity and storage conditions of the RNA template. Make new RNA dilutions from the stock if necessary.</li> </ul>
RNA template contains strong secondary structures	<ul style="list-style-type: none"> <li>Perform a predenaturation step on the template before cDNA synthesis. See Section 4.3.</li> </ul>
Incorrect temperature in cDNA synthesis reaction	<ul style="list-style-type: none"> <li>The recommended temperature in cDNA synthesis step is 37°C. It can be optimized between 37–48°C if necessary.</li> </ul>
RT-PCR product too long	<ul style="list-style-type: none"> <li>The length of the amplicon should be between 50 and 250 bp. The optimal length is 100–150 bp.</li> </ul>
Insufficient activation of the hot start <i>Tbr</i> DNA polymerase	<ul style="list-style-type: none"> <li>Make sure 95°C 15 min was used for the initial denaturation step.</li> <li>Make sure the cycler block temperature is accurate.</li> </ul>
Insufficient extension time for the amplicon size	<ul style="list-style-type: none"> <li>Increase extension time.</li> </ul>
qPCR primer or probe design not optimal	<ul style="list-style-type: none"> <li>Check primer or probe design. See Section 5.1.</li> </ul>
qPCR primer or probe concentration too low	<ul style="list-style-type: none"> <li>Increase qPCR primer concentration (to a maximum of 1 µM each). 0.25 µM probe concentration is usually sufficient.</li> </ul>

Annealing temperature too high in qPCR	<ul style="list-style-type: none"> <li>• Use a gradient to optimize the annealing temperature.</li> <li>• Decrease the annealing temperature in 2°C decrements if no gradient feature is available.</li> </ul>
PCR protocol not optimal	<ul style="list-style-type: none"> <li>• Make sure you are using the recommended PCR protocol. If necessary, optimize using the recommended protocol as a starting point.</li> </ul>
<b>Increase in fluorescence signal in negative (no RT) control</b>	
Contaminating genomic DNA in RNA preparation	<ul style="list-style-type: none"> <li>• Treat the starting RNA template with DNase I before cDNA synthesis.</li> <li>• Redesign qPCR primers or probe to prevent amplification of genomic DNA. See Section 5.1.</li> </ul>
PCR products from a previous run contaminating the reaction	<ul style="list-style-type: none"> <li>• Perform UNG treatment before PCR cycling.</li> </ul>
Weak RT activity of the DNA polymerase	<ul style="list-style-type: none"> <li>• Due to a weak RT activity of most DNA polymerases, short RNA target sequences might be reverse transcribed during qPCR step.</li> <li>• If contamination and amplification of genomic DNA can be excluded, you can ignore this no RT control.</li> </ul>
<b>Normal fluorescence signal, but low efficiency</b>	
Missing components or pipetting error	<ul style="list-style-type: none"> <li>• Check the assembly of the reactions.</li> <li>• Check the concentrations and storage conditions of the reagents.</li> </ul>
Primer and probe design not optimal or low template concentration	<ul style="list-style-type: none"> <li>• Check primer and probe design and template stock concentration.</li> <li>• Increase template amount (up to 10 % of the qPCR reaction volume).</li> </ul>
Inhibitors from the sample affecting reaction	<ul style="list-style-type: none"> <li>• Repurify RNA.</li> <li>• Reduce template amount. The volume of the cDNA template should not exceed 10 % of the qPCR reaction volume.</li> </ul>
<b>Non-linear correlation between Ct and log of template amount in the standard curve</b>	
RNA template dilution inaccurate	<ul style="list-style-type: none"> <li>• Remake dilution series and make sure the samples are well mixed.</li> </ul>
RNA template amount too high	<ul style="list-style-type: none"> <li>• Do not exceed 1 µg of starting RNA template.</li> </ul>
RNA template amount too low	<ul style="list-style-type: none"> <li>• Increase the amount of RNA template (max 1 µg).</li> </ul>
RT-PCR product too long	<ul style="list-style-type: none"> <li>• The length of the amplicon should be between 50 and 250 bp. The optimal length is 100–150 bp.</li> </ul>
cDNA template volume too high	<ul style="list-style-type: none"> <li>• Reduce template amount. The volume of the cDNA template should not exceed 10 % of the qPCR reaction volume.</li> <li>• Increase qPCR reaction volume.</li> </ul>
cDNA template volume too low	<ul style="list-style-type: none"> <li>• Increase template amount (up to 10 % of qPCR reaction volume).</li> </ul>
Insufficient activation of the hot-start DNA polymerase	<ul style="list-style-type: none"> <li>• Make sure 95°C 15 min was used for the initial reactivation/denaturation step in qPCR.</li> <li>• Make sure cycler block temperature is accurate.</li> </ul>
Insufficient denaturation of the template	<ul style="list-style-type: none"> <li>• Make sure 95°C 15 min is used for the initial denaturation step</li> <li>• Make sure cycler block temperature is accurate.</li> </ul>
Co-amplification of primer-dimers with the specific product	<ul style="list-style-type: none"> <li>• Redesign primers.</li> </ul>

qPCR primer/probe design or concentration not optimal	<ul style="list-style-type: none"> <li>Check primer/probe design and concentration. See Section 5.1.</li> </ul>
<b>Low signal when using ROX normalization</b>	
High ROX passive reference fluorescence intensity	<ul style="list-style-type: none"> <li>Use lower ROX concentration. See recommended concentrations in Table 2.</li> </ul>
<b>High signal when using ROX normalization</b>	
Low ROX passive reference fluorescence intensity	<ul style="list-style-type: none"> <li>Use higher ROX concentration. See recommended concentrations in Table 5.</li> </ul>
Yellow dye in the sample buffer decreases ROX intensity	<ul style="list-style-type: none"> <li>Use higher ROX concentration. See recommended concentrations in Table 5.</li> </ul>
<b>Abnormal appearance of amplification curves when ROX normalization is used</b>	
Color calibration not accurate. Fluorescence intensity from one channel affects intensity in another channel	<ul style="list-style-type: none"> <li>Verify color calibration according to instrument instructions.</li> </ul>

## Appendix I: general molecular biology data

Table 10. Spectrophotometric conversions for nucleic acid templates.

<b>1 A<sub>260</sub> unit*</b>	<b>Concentration (µg/ml)</b>
Double-stranded DNA	50
Single-stranded DNA	33
Single-stranded RNA	40

\* Absorbance at 260 nm = 1 (1 cm detection path).

Table 11. Molar conversions for nucleic acid templates.

<b>Nucleic acid</b>	<b>Size</b>	<b>pmol/µg</b>	<b>Copies/µg*</b>
1 kb DNA	1 000 bp	1.52	$9.1 \times 10^{11}$
pUC19DNA	2 686 bp	0.57	$3.4 \times 10^{11}$
Lambda DNA	48 502 bp	0.03	$1.8 \times 10^{10}$
<i>Escherichia coli</i>	$4.7 \times 10^6$ bp	$3.2 \times 10^{-4}$	$1.9 \times 10^9$
Human	$3.2 \times 10^9$ bp	$4.7 \times 10^{-7}$	$2.8 \times 10^5$

\* For single-copy genes.

## Product use limitation

This product has been developed and is sold exclusively for research purposes and in vitro use only. This product has not been tested for use in diagnostics or drug development, nor is it suitable for administration to humans or animals.

## Trademark and patent notices; label licenses

DyNAmo is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

ROX is a trademark of Life Technologies Corporation and its affiliated companies.

Applied Biosystems is a registered trademark of Life Technologies Corporation and its affiliated companies.

SYBR is a registered trademark of Life Technologies Corporation and its affiliated companies.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.

Use of this product is covered by US Patent No. 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California

Use of this product in a passive reference method is covered by the following U.S. Patent: 5,928,907 (claim numbers 12-24, 27-28) and corresponding patent claims outside the US. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product in a passive reference method for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

The quality system of Finnzymes Oy is certified according to standard SFS-EN ISO9001:2008.

**Thermo Scientific DyNAmo Probe 2-Step qRT-PCR Kit  
Technical Manual**



© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

v3\_10.2011

Technical support:  
US: [techservice.genomics@thermofisher.com](mailto:techservice.genomics@thermofisher.com)

Europe, Asia, Rest of World:  
[techservice.emea.genomics@thermofisher.com](mailto:techservice.emea.genomics@thermofisher.com)

**Thermo**  
SCIENTIFIC