

Thermo Scientific Phusion High-Fidelity PCR Kit

F-553S/L, 50 U/200 U

Store at -20°C

Thermo
SCIENTIFIC

1. Introduction

Thermo Scientific Phusion High-Fidelity DNA Polymerase offers extreme performance for all major PCR applications. Incorporating an exciting technology, Phusion® DNA Polymerase brings together a novel *Pyrococcus*-like enzyme with a processivity-enhancing domain. Phusion DNA Polymerase generates long templates with an accuracy and speed previously unattainable with a single enzyme, even on the most difficult templates. The extreme fidelity makes Phusion DNA Polymerase a superior choice for cloning. The error rate of Phusion DNA Polymerase in Phusion HF Buffer is 4.4×10^{-7} when determined with a modified *lacI*-based method¹. It is approximately 50-fold lower than that of *Thermus aquaticus* DNA polymerase and 6-fold lower than that of *Pyrococcus furiosus* DNA polymerase.

The Phusion High-Fidelity PCR Kit includes lambda DNA control template and primers for 1.3 kb and 10 kb amplicons. The template amount is sufficient for performing 20 control reactions in 50 µl volume or 50 control reactions in 20 µl volume.

IMPORTANT NOTES

- Use 98°C for denaturation (see 5.1 & 5.2).
- The annealing rules are different from many common DNA polymerases (such as *Taq* DNA polymerases). Read Section 5.3 carefully.
- Use 15–30 s/kb for extension. Do not exceed 1 min/kb (see 5.4).
- Use Phusion DNA Polymerase at 0.5–1.0 U per 50 µl reaction volume. Do not exceed 2 U/50 µl (see 4.1).
- Use 200 µM of each dNTP. Do not use dUTP (see 4.3).
- Note: Phusion DNA Polymerases produce blunt end DNA products.

2. Kit components

Component	Concentration	F-553S	F-553L
Phusion DNA Polymerase	2 U/µl	50 U	200 U
5x Phusion HF buffer		1 x 1.5 ml	3 x 1.5 ml
5x Phusion GC buffer		1 x 1.5 ml	1 x 1.5 ml
dNTP mix	10 mM each	100 µl	400 µl
MgCl ₂ solution	50 mM	1.5 ml	1.5 ml
Control lambda template	0.5 ng/µl	40 µl	40 µl
1.3 kb primers	4 µM each	50 µl	50 µl
10 kb primers	4 µM each	50 µl	50 µl
DNA size standard		200 µl	400 µl
DMSO		0.5 ml	0.5 ml

Material safety data sheet (MSDS) is available at www.thermoscientific.com/fzmsds.

3. Guidelines for Phusion DNA Polymerase

Phusion DNA Polymerase (2U/µl) is provided with 5x Phusion HF Buffer and 5x Phusion GC Buffer. Both buffers contain 1.5 mM MgCl₂ at final reaction concentrations. Separate tubes of DMSO and 50 mM MgCl₂ solutions are provided for further optimization.

3.1 Basic reaction conditions for DNA amplification

Carefully mix and centrifuge all tubes before opening to ensure homogeneity and improve recovery. PCR reactions should be set up on ice. Phusion DNA Polymerase should be pipetted carefully and gently as the high glycerol content (50 %) in the storage buffer may otherwise lead to pipetting errors. It is critical that Phusion DNA Polymerase is the last component added to the PCR mixture, since the enzyme exhibits 3'→5' exonuclease activity that can degrade primers in the absence of dNTPs.

Due to the nature of Phusion DNA Polymerase, the optimal reaction conditions may differ from PCR protocols for standard DNA polymerases. Due to the high salt concentration in the reaction buffer, Phusion DNA Polymerase tends to work better at elevated denaturation and annealing temperatures. Please pay special attention to the conditions listed in section 5 when running your reactions. Following the guidelines will ensure optimal enzyme performance.

Table 1. Pipetting instructions (add items in this order).

Component	50 µl react.	20 µl react.	Final conc.
H ₂ O	add to 50 µl	add to 20 µl	
5x Phusion HF Buffer*	10 µl	4 µl	1x
10 mM dNTPs	1 µl	0.4 µl	200 µM each
primer A**	x µl	x µl	0.5 µM
primer B**	x µl	x µl	0.5 µM
template DNA	x µl	x µl	
(DMSO***, optional)	(1.5 µl)	(0.6 µl)	(3 %)
Phusion DNA Polymerase	0.5 µl	0.2 µl	0.02U/µl

* Optionally 5x Phusion GC Buffer can be used, see section 4.2 for details.

** The recommendation for final primer concentration is 0.5 µM, but it can be varied in the range of 0.2–1.0 µM, if needed.

*** Addition of DMSO is recommended for GC-rich amplicons. DMSO is not recommended for amplicons with very low GC % or amplicons that are >20 kb.

Table 2. Cycling instructions.

Cycle step	2-step protocol		3-step protocol		Cycles
	Temp.	Time	Temp.	Time	
Initial denaturation	98°C	30 s	98°C	30 s	1
Denaturation Annealing (see 5.3) Extension (see 5.4)	98°C	5–10 s	98°C	5–10 s	25–35
	–	–	X°C	10–30 s	
	72°C	15–30 s/1 kb	72°C	15–30 s/1 kb	
Final extension	72°C	5–10 min hold	72°C	5–10 min hold	1
	4°C		4°C		

4. Notes about reaction components

4.1 Enzyme

The optimal amount of enzyme depends on the amount of template and the length of the PCR product. Usually 1 unit of Phusion DNA Polymerase per 50 µl reaction volume gives good results, but optimal amounts can range from 0.5 to 2 units per 50 µl reaction depending on amplicon length and difficulty. **Do not exceed 2 U/50 µl (0.04 U/µl), especially for amplicons that are > 5kb.**

When cloning fragments amplified with Phusion DNA Polymerase, blunt end cloning is recommended. If TA cloning is required, it can be performed by adding A overhangs to the blunt PCR product with Thermo Scientific *Taq* DNA Polymerase, for example. However, before adding the overhangs it is very important to remove all the Phusion DNA Polymerase by purifying the PCR product carefully. Any remaining Phusion DNA Polymerase will degrade the A overhangs, creating blunt ends again. A detailed protocol for TA cloning of DNA fragments amplified with any of the Phusion DNA polymerases can be found on website: www.thermoscientific.com/pcrcloning.

4.2 Buffers

Two buffers are provided with the enzyme: 5x Phusion HF Buffer (F-518) and 5x Phusion GC Buffer (F-519). The error rate of Phusion DNA Polymerase in HF Buffer (4.4×10^{-7}) is lower than that in GC Buffer (9.5×10^{-7}). Therefore, the HF Buffer should be used as the default buffer for high-fidelity amplification. However, GC Buffer can improve the performance of Phusion DNA Polymerase on some difficult or long templates, such as GC-rich templates or those with complex secondary structures. For applications such as microarray or DHPLC, where the DNA templates need to be free of detergents, detergent-free reaction buffers (F-520, F-521) are available for Phusion DNA Polymerases.

4.3 Mg²⁺ and dNTP

The concentration of Mg²⁺ is critical since Phusion DNA Polymerase is a magnesium dependent enzyme. Excessive Mg²⁺ stabilizes the DNA double strand and prevents complete denaturation of DNA. Excess Mg²⁺ can also stabilize spurious annealing of primers to incorrect template sites and decrease specificity. Conversely, inadequate Mg²⁺ may lead to lower product yield. The optimal Mg²⁺ concentration also depends on the dNTP concentration, the specific template DNA and the sample buffer composition. In general, the optimal Mg²⁺ concentration is 0.5 to 1 mM over the total dNTP concentration for standard PCR. If the primers and/or template contain chelators such as EDTA or EGTA, the apparent Mg²⁺ optimum may be shifted to higher concentrations. If further optimization is needed, increase Mg²⁺ concentration in 0.2 mM steps.

High quality dNTPs should be used for optimal performance with Phusion DNA Polymerase. The polymerase cannot read uracil derivatives or inosine in the template strand so the use of these analogues or primers containing them is not recommended. Due to the high processivity of Phusion DNA Polymerase there is no advantage in increasing dNTP concentrations. For optimal results always use 200 µM of each dNTP.

4.4 Template

General guidelines for low complexity DNA (e.g. plasmid, lambda or BAC DNA) are: 1 pg–10 ng per 50 µl reaction volume. For high complexity genomic DNA, the amount of DNA template should be 50–250 ng per 50 µl reaction volume. If cDNA synthesis reaction mixture is used as a source of template, the volume of the template should not exceed 10 % of the final PCR reaction volume.

4.5 PCR additives

The recommended reaction conditions for GC-rich templates include 3 % DMSO as a PCR additive, which aids in the denaturing of templates with high GC contents. For further optimization DMSO should be varied in 2 % increments. In some cases DMSO may also be required for supercoiled plasmids to relax for denaturation. Other PCR additives such as formamide, glycerol, and betaine are also compatible with Phusion DNA Polymerase.

If high DMSO concentration is used, the annealing temperature must be decreased, as DMSO affects the melting point of the primers. It has been reported that 10 % DMSO decreases the annealing temperature by 5.5–6.0°C².

5. Notes about cycling conditions

5.1 Initial denaturation

Denaturation should be performed at 98°C. Due to the high thermostability of Phusion DNA Polymerase even higher than 98°C denaturation temperatures can be used. We recommend a 30-second initial denaturation at 98°C for most templates. Some templates may require longer initial denaturation time and the length of the initial denaturation time can be extended up to 3 minutes.

5.2 Denaturation

Keep the denaturation time as short as possible. Usually 5–10 seconds at 98°C is enough for most templates. **Note:** The denaturation time and temperature may vary depending on the ramp rate and temperature control mode of the cyclor.

5.3 Primer annealing

The optimal annealing temperature for Phusion Hot Start DNA Polymerase may differ significantly from that of *Taq*-based polymerases. Always use the T_m calculator and instructions on website: www.thermoscientific.com/pcrwebtools to determine the T_m values of primers and optimal annealing temperature. As a basic rule, for primers > 20 nt, anneal for 10–30 seconds at a T_m+3°C of the lower T_m primer. For primers ≤ 20 nt, use an annealing temperature equal to the T_m of the lower T_m primer. If necessary, use a temperature gradient to find the optimal annealing temperature for each template-primer pair combination. The annealing gradient should extend up to the extension temperature (two-step PCR). A 2-step protocol is recommended when primer T_m values are at least 69°C (> 20 nt) or 72°C (≤ 20 nt) when calculated with Thermo Scientific's T_m calculator. In the 2-step protocol the combined annealing/extension step should be performed at 72°C even when the primer T_m is > 72°C.

5.4 Extension

The extension should be performed at 72°C. Extension time depends on amplicon length and complexity. For low complexity DNA (e.g. plasmid, lambda or BAC DNA) use an extension time of 15 seconds per 1 kb. For high complexity genomic DNA 30 seconds per 1 kb is recommended. For some cDNA templates, the extension time can be increased up to 40 seconds per 1 kb to obtain optimal results.

6. Amplifying control template

6.1 Reaction conditions

Table 3. Pipetting instructions for control reactions.

Component	50 µl react.	20 µl react.	Final conc.
H ₂ O	34 µl	13.6 µl	
5x Phusion HF Buffer	10 µl	4 µl	1x
10 mM dNTPs	1 µl	0.4 µl	200 µM each
Primers*	2.5 µl	1 µl	0.2 µM
Control template DNA	2 µl	0.8 µl	
Phusion DNA Polymerase	0.5 µl	0.2 µl**	0.02U/µl

* Either the 1.3 kb primer set or 10 kb primer set.

** Dilution of polymerase should be made to 1x reaction buffer to avoid pipetting errors.

6.2 Cycling conditions

A separate cycling protocol is given for both 1.3 kb and 10 kb control amplicons. Alternatively, both control reactions can be amplified simultaneously using the 10 kb cycling protocol.

Table 4. Cycling conditions for 1.3 kb control fragment (2-step protocol).


Cycle step	Temp.	Time	Cycles
Initial denaturation	98°C	1 min	1
Denaturation Annealing / Extension	98°C	5 s	25
	72°C	20 s	
Final extension	72°C	10 min hold	1
	10°C		

Table 5. Cycling conditions for 10 kb fragment (3-step protocol). This program can also be used if both control reactions are amplified simultaneously.

Cycle step	Temp.	Time	Cycles
Initial denaturation	98°C	1 min	1
Denaturation Annealing Extension	98°C	5 s	25
	60°C	15 s	
	72°C	2 min 30 s	
Final extension	72°C	10 min hold	1
	10°C		

The cycling protocols above are recommendations. If you wish to run these controls together or with your experimental samples, please note that the controls have been shown to work in a variety of conditions. The 1.3 kb control has been successfully amplified with both 2- and 3-step protocols with extension times ranging from 15 s to 5 min, and cycle numbers ranging from 20 to 30. The 10 kb control has been successfully amplified with 3-step protocol with extension times ranging from 2 min to 5 min, and cycle numbers ranging from 20 to 30.

6.3 Analysis of the control reactions



In the image on the left both control reactions have been run on an ethidium bromide stained agarose gel (1% SeaKem LE agarose in TAE buffer). For this run 15 µl of loading dye was added to the 50 µl control PCR reactions, and 5 µl of the resulting mixtures were loaded on the gel.

Lane 1. DNA size standard
Lane 2. 1.3 kb control amplicon
Lane 3. 10 kb control amplicon
Lane 4. DNA size standard

After running your control reactions on a gel, compare the results to the image on the left to check for specificity and efficiency of the reactions.

7. Troubleshooting

No product at all or low yield	
<ul style="list-style-type: none"> Repeat and make sure that there are no pipetting errors. Use fresh high quality dNTPs. Do not use dNTP mix that contain dUTP or dITP or primers that contain uracil or inosine. Titrate template amount. Template DNA may be damaged. Use carefully purified template. Increase extension time. Increase the number of cycles. Optimize annealing temperature. Optimize enzyme concentration. 	<ul style="list-style-type: none"> Titrate DMSO (2–8 %) in the reaction (see section 4.5). Denaturation temperature may be too low. Optimal denaturation temperature for most templates is 98°C or higher. Optimize denaturation time. Check the purity and concentration of the primers. Check primer design. Try using the alternative GC buffer (see section 4.2). If DNA is not carefully purified, inhibitors may be present – decrease the amount of DNA.
Non-specific products - High molecular weight smears	
<ul style="list-style-type: none"> Decrease enzyme concentration (see section 4.1). Make sure the extension time used was not too long (see section 5.4). Titrate template amount. Reduce the total number of cycles. 	<ul style="list-style-type: none"> Increase annealing temperature or try 2-step protocol (see section 5.3). Vary denaturation temperature (see section 5.2). Optimize Mg²⁺-concentration. Decrease primer concentration.
Non-specific products - Low molecular weight discrete bands	
<ul style="list-style-type: none"> Increase annealing temperature (see section 5.3). Shorten extension time (see section 5.4). Decrease enzyme concentration. 	<ul style="list-style-type: none"> Optimize Mg²⁺-concentration. Titrate template amount. Decrease primer concentration. Design new primers.

8. Component specifications

8.1 Phusion® High-Fidelity DNA Polymerase (F-530)

Thermostable Phusion DNA Polymerase is purified from an *E. coli* strain expressing the cloned Phusion DNA Polymerase gene. Phusion DNA Polymerase possesses the following activities: 5'→3' DNA polymerase activity and 3'→5' exonuclease activity. Phusion DNA Polymerase is free of contaminating endo- and exonucleases.

Storage buffer: 20 mM Tris-HCl (pH 7.4 at 25°C), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, 200 µg/ml BSA and 50 % glycerol. **Unit definition:** One unit is defined as the amount of enzyme that will incorporate 10 nmoles of dNTPs into acid-insoluble form at 74°C in 30 minutes under the stated assay conditions.

Unit assay conditions: Incubation buffer: 25 mM TAPS-HCl, pH 9.3 (at 25°C), 50 mM KCl, 2 mM MgCl₂, 1 mM β-mercaptoethanol, 100 µM dCTP, 200 µM each dATP, dGTP, dTTP. Incubation procedure: 20 µg activated calf thymus DNA and 0.5 µCi [α -³²P] dCTP are incubated with 0.1 units of DNA polymerase in 50 µl incubation buffer at 74°C for 10 minutes. The amount of incorporated dNTPs is determined by trichloroacetic acid precipitation.

DNA amplification assay: Performance in PCR is tested by the amplification of a 7.5 kb fragment of genomic DNA and a 20 kb fragment of lambda DNA.

Exonuclease contamination assay: Incubation of 10 U for 4 hours at 72°C in 50 µl assay buffer with 1 µg sonicated [³H] ssDNA (2x10⁵ cpm/µg) released < 1 % of radioactivity.

Endonuclease contamination assay: No endonuclease activity is observed after incubation of 10 U of DNA polymerase with 1 µg of λ DNA in assay buffer at 72°C for 4 hours.

8.2 5x Phusion® HF Buffer (F-518)

The 5x Phusion HF Buffer contains 7.5 mM MgCl₂, which provides 1.5 mM MgCl₂ in final reaction conditions.

8.3 5x Phusion® GC Buffer (F-519)

The 5x Phusion GC Buffer contains 7.5 mM MgCl₂, which provides 1.5 mM MgCl₂ in final reaction conditions.

8.4 dNTP mix (F-560)

The dNTP mix is a premixed ready-to-use solution consisting of the following compounds: dATP, dGTP, dCTP and dTTP dissolved in H₂O at 10 mM each.

8.5 50 mM MgCl₂ solution (F-510MG)

Both Phusion Buffers supply 1.5 mM MgCl₂ at final reaction conditions. If higher MgCl₂ concentrations are desired, use 50 mM MgCl₂ solution to increase the MgCl₂ titer. Using the following equation, you can calculate the volume of 50 mM MgCl₂ needed to attain the final MgCl₂ concentration:

[desired mM Mg] - [1.5 mM] = µl to add to a 50 µl reaction.

For example, to increase the MgCl₂ concentration to 2.0 mM, add 0.5 µl of the 50 mM MgCl₂ solution. Because the PCR reactions can be quite sensitive to changes in the MgCl₂ concentration, it is recommended that the 50 mM MgCl₂ stock solution is diluted 1:5 (to 10 mM) to minimize pipetting errors.

8.6 Lambda DNA control template (F-304K)

The control template is bacteriophage lambda DNA (GenBank access number J02459, 48 502 bp). The concentration is 0.5 ng/µl in TE buffer.

8.7 1.3 kb control primer mix (F-535)

This component is a mix of primers for amplification of a 1.3 kb fragment of lambda DNA. Each primer concentration is 4 µM in H₂O.

Primer #1 (27-mer)
5'-GTC ACC AGT GCA GTG CTT GAT AAC AGG-3'
Melting point: 71.0°C
Coordinates in lambda DNA: 30 006–30 032

Primer #2 (28-mer)
5'-GAT GAC GCA TCC TCA CGA TAA TAT CCG G-3'
Melting point: 73.2°C
Coordinates in lambda DNA: 31 325–31 352

8.8 10 kb control primer mix (F-536)

This component is a mix of primers for amplification of a 10 kb fragment of lambda DNA. Each primer concentration is 4 μ M in H₂O.

Primer #1 (22-mer)
5'-CAG TGC AGT GCT TGA TAA CAG G-3'
Melting point: 63.5°C
Coordinates in lambda DNA: 30 011–30 032

Primer #2 (20-mer)
5'-GTA GTG CGC GTT TGA TTT CC-3'
Melting point: 63.3°C
Coordinates in lambda DNA: 40 024–40 043

8.9 Ready-to-use DNA standard (F-303SD)

This DNA standard is a mix of lambda DNA HindIII digest and bacteriophage ϕ X174 DNA HaeIII digest, each at 50 ng/ μ l (100 ng/ μ l total). It is supplied in 8 mM Tris-HCl (pH 8.0), 12 mM EDTA, 12 % glycerol and 0.012 % (w/v) bromophenol blue dye.

The DNA standard solution contains 19 fragments of the following sizes and mass amounts (per 10 μ l):



Fragment	Base pairs	DNA amount ng/10 μ l
1	23 130	238
2	9 416	97
3	6 557	68
4	4 361	45
5	2 322	24
6	2 027	21
7	1 353	126
8	1 078	100
9	872	81
10	603	56
11	564*	6
12	310	29
13a	281	26
13b	271	25
14	234	22
15	194	18
16a	125*	1
16b	118	11
17	72	7

Note: The cohesive areas of fragments 1 and 4 can be separated by heating at 65°C for 5 minutes. For daily use the marker can be stored at +4°C (at least one month). The marker is stable at -20°C for at least one year.

* Due to the low amount of DNA these bands are almost invisible.

8.10 Dimethyl sulfoxide DMSO, 100 % (F-515)

Note: The freezing point of DMSO is 18–19°C, so it does not melt on ice.

9. References

1. Frey M. & Suppmann B. (1995) *Biochemica* 2: 34–35.
2. Chester N. & Marshak D.R. (1993) *Analytical Biochemistry* 209: 284–290.

Shipping and storage

Phusion DNA Polymerase is shipped on gel ice. Upon arrival, store the components at -20°C.

Technical support

US: techservice.genomics@thermofisher.com
Europe, Asia, Rest of World:
techservice.emea.genomics@thermofisher.com

Web: www.thermoscientific.com/phusion
Tm-calculator: www.thermoscientific.com/pcrwebtools

Product use limitation

This product has been developed and is sold exclusively for research purposes and in vitro use only. This product has not been tested for use in diagnostics or drug development, nor are they suitable for administration to humans or animals.

Notice to Purchaser: Limited License

The purchase price of this product includes a limited, non-transferable license under national patents from EP 0 547 920 B1, owned by New England Biolabs, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

The purchase price of this product includes a limited, nontransferable license under U.S. and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

Designed and manufactured according to certified ISO9001:2008 processes.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

v1_07.2011